

Allyl D-Glucuronate

Art. ID TRC-A555400-500MG
Unit 500 mg
Deliverydetails Dry ice shipment

Description

Category: Carbohydrates and Derivatives, Glucuronides /// Appearance: White to Off-White Solid /// Application Notes: It is used in efficient synthesis of 1beta-O-acyl glucuronides via regioselective acylation of allyl or benzyl D-glucuronate. alphasyl glucuronides are known metabolites of important drugs. /// References: Kaspersen, F., et al.: Xenobiotica, 17, 1451

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Allyl D-Glucuronate	[188717-04-6]				