

### **3-Aminopropyltriethoxy?silane**

Art. ID TRC-A628785-100G  
Unit 100 g  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2735 II E2

#### Description

Appearance: Clear Colourless Oil /// Application Notes: 3-Aminopropyltriethoxy?silane is a reagent in the preparation of ultrastable organic dye nanoparticles in living cells. Also used in the preparation of monoclonal antibody-functionalized mesoporous silica nanoparticles packaged with Mifepristone (M343975), a progesterone receptor antagonist with partial agonist activity. Abortifacient. /// References: Xu, R. et al.: *Biomat.*, 93, 38 (2016)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silquest A-1100	[919-30-2]				