

(4-Acetyl-L-Serine)Leuprorelin TFA Salt

Art. ID TRC-A187985-25MG
Unit 25 mg
Deliverydetails No Dangerous Good

Description

Appearance: White Solid /// Application Notes: (4-Acetyl-L-Serine)Leuprorelin TFA Salt is an impurity of Leuprolide (L330600), which is a synthetic nonapeptide agonist analog of LH-RH. /// References: Vilchez-Martinez, J.A., et al.: Biochem. Biophys. Res. Commun., 59, 1226 (1974)| Sennello, L.T., et al.: J. Pharm. Sci., et al.: 75, 158 (1986)| Wheeler, J.M., et al.: Am.J. Obstet. Gynecol., 167, 1367 (1992)| Eri, L.M., et al.: J. Urol., 150, 359 (1993)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(4-Acetyl-L-Serine)Leup rorelin TFA Salt					