

D-Alanine Isopropyl Ester Hydrochloride

Art. ID TRC-A510820-50G
Unit 50 g
Deliverydetails No Dangerous Good

Description

Appearance: White to Off-White Solid /// Application Notes: D-Alanine Isopropyl Ester Hydrochloride is a chiral reagent used in the synthesis of nucleotide prodrugs such as (PSI-7977) for the treatment of hepatitis C. Also used in the preparation of synthetic sweeteners. /// References: Sofia, M. et al.: J. Med. Chem., 53, 7202 (2010)| Mazur, R. et al.: J. Med. Chem., 16, 1284 (1973)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	D-Alanine Isopropyl Ester Hydrochloride	[39613-92-8]				