

A77 1726

Art. ID TRC-A771726-500MG
Unit 500 mg
Deliverydetails De Minimis UN2811 III E1

Description

Category: Metabolites and Impurities, Tyrosine Kinase Inhibitors, Pharmaceuticals, Intermediates and Fine Chemicals /// Appearance: White Solid /// Application Notes: A77 1726 (Leflunomide EP Impurity B) is the active metabolite of Leflunomide, a potent disease-modifying antirheumatic drug used in the treatment of rheumatoid arthritis. LEF-M interferes with dendritic cell function. The second Phase III study evaluating the efficacy of Genzyme's teriflunomide therapy against relapsing forms of multiple sclerosis found the once-daily oral drug to be no less effective than injectable interferon beta 1a (Rebif®). /// References: Kirsch, B.M., et al.: Arthritis Research and Therapy, 7, 3, R694, (2005)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	A77 1726	[163451-81-8]				