

**Lead-Tin Solder, Cast, Certified Reference Material (CRM), disc 40 mm x 15 mm**

Art. ID MBH-91XS30PR2C  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,06	%		
	Aluminium (Al)	[7429-90-5]	<0.0005	%		
	Arsenic (As)	[7440-38-2]	0,028	%		
	Gold (Au)	[7440-57-5]	0,0017	%		
	Bismuth (Bi)	[7440-69-9]	0,158	%		
	Cadmium (Cd)	[7440-43-9]	0,006	%		
	Copper (Cu)	[7440-50-8]	0,095	%		
	Iron (Fe)	[7439-89-6]	0,009	%		
	Nickel (Ni)	[7440-02-0]	0,0077	%		
	Antimony (Sb)	[7440-36-0]	0,619	%		
	Tin (Sn)	[7440-31-5]	30,17	%		
	Zinc (Zn)	[7440-66-6]	0,016	%		