

**Lead-Tin Solder, Cast, Certified Reference Material (CRM), disc 40 mm x 15 mm**

Art. ID MBH-91XS63Bi1A  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,0592	%		
	Arsenic (As)	[7440-38-2]	(<0.002)	%		
	Gold (Au)	[7440-57-5]	0,074	%		
	Bismuth (Bi)	[7440-69-9]	0,597	%		
	Cadmium (Cd)	[7440-43-9]	0,0095	%		
	Copper (Cu)	[7440-50-8]	0,105	%		
	Iron (Fe)	[7439-89-6]	0,0204	%		
	Indium (In)	[7440-74-6]	0,0067	%		
	Nickel (Ni)	[7440-02-0]	0,0131	%		
	Antimony (Sb)	[7440-36-0]	0,47	%		
	Tin (Sn)	[7440-31-5]	61,9	%		
	Tellurium (Te)	[13494-80-9]	0,0012	%		