

Phosphatidylcholine, PC (egg), Purity: >98%

Art. ID LA-37-0100-9
Unit 100 mg
Deliverydetails No Dangerous Good /not restricted

Description

Synonyms: PC (egg) / L-?-phosphatidylcholine (egg)

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|----------------------------------|------------|---------------------|------|--------|--------|
| | Phosphatidylcholine, PC (egg) | | | | | |