

p-Aminophenyl Arsenoxide

Art. ID TRC-A622500-100MG
Unit 100 mg
Deliverydetails De Minimis

Description

Category: Aromatics, Miscellaneous Reagents /// Appearance: Off-White Solid /// Application Notes: p-Aminophenyl Arsenoxide is a useful reagent for the study of 2-oxoacid dehydrogenase multienzyme complexes. /// References: Stevenson, et al.: Biochemistry, 17, 11, 2191 (1978)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	p-Aminophenyl Arsenoxid	[1122-90-3]				
	e					