

Phosphatidylinositol, PI (plant) NH₄⁺ salt, Purity: >98%, solution in Chloroform, concentration: Varies

Art. ID LA-37-0130-7
Unit 25 mg
Deliverydetails De Minimis UN1888 III E1

Description

Synonyms: PI (plant) Ammonium salt / L-?-phosphatidylinositol (soybean)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphatidylinositol, P		Varies			
	I (plant) NH ₄ ⁺ salt					