

**1,2-Dihexanoyl-sn-Glycero-3-Phosphatidylcholine, PC(6:0/6:0), Purity: >99%, solution in Chloroform, concentration: Varies**

Art. ID LA-37-0600-9  
Unit 100 mg  
Deliverydetails De Minimis UN1888 III E1

Description

Synonyms: 6:0 PC / 1,2-dicaproyl-sn-glycero-3-phosphocholine / [(2R)-2,3-di(hexanoyloxy)propyl] 2-trimethylazaniumylethyl phosphate / / [(2R)-2,3-bis(1-oxohexoxy)propyl] 2-trimethylammonioethyl phosphate

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2-Dihexanoyl-sn-Glycero-3-Phosphatidylcholine	[53892-41-4]	Varies			