

### **Custom Solution - Spiking Standard 2R - 4 components in 2% nitric acid**

Art. ID LM24-CUS-70181  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

#### Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Magnesium (Mg)	[7439-95-4]	10000	mg/L		
	Potassium (K)	[7440-09-7]	10000	mg/L		
	Sodium (Na)	[7440-23-5]	10000	mg/L		
	Calcium (Ca)	[7440-70-2]	10000	mg/L		