

## **Custom Solution - 8 components in 5% HNO<sub>3</sub>, trace HF**

Art. ID LM24-CUS-70184

Unit 50 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

### **Description**

Conservation 24 months (12 months after opening)

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Arsenic (As)      | [7440-38-2] | 2000                | mg/L |        |        |
|                  | Cadmium (Cd)      | [7440-43-9] | 50                  | mg/L |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | 1000                | mg/L |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 1000                | mg/L |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 1000                | mg/L |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 1000                | mg/L |        |        |
|                  | Tin (Sn)          | [7440-31-5] | 1000                | mg/L |        |        |
|                  | Zinc (Zn)         | [7440-66-6] | 1000                | mg/L |        |        |