

Custom Solution - 5 components in 5% HNO₃

Art. ID LM24-CUS-70185

Unit 50 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	50	mg/L		
	Cadmium (Cd)	[7440-43-9]	100	mg/L		
	Copper (Cu)	[7440-50-8]	2000	mg/L		
	Lead (Pb)	[7439-92-1]	2000	mg/L		
	Zinc (Zn)	[7440-66-6]	5000	mg/L		