

Aluminium base, CRM, Alloy: 3003H, disc 50x30mm

Art. ID DHJ-58A-ZB23007
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Copper (Cu) | [7440-50-8] | 0,109 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,38 | % | | |
| | Magnesium (Mg) | [7439-95-4] | 0,018 | % | | |
| | Manganese (Mn) | [7439-96-5] | 1,07 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,0036 | % | | |
| | Lead (Pb) | [7439-92-1] | <0,0010 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,128 | % | | |
| | Tin (Sn) | [7440-31-5] | <0,0010 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,016 | % | | |
| | Vanadium (V) | [7440-62-2] | 0,014 | % | | |
| | Zinc (Zn) | [7440-66-6] | 0,004 | % | | |