

Aluminium base, RM, Alloy: A380, disc 45x35mm

Art. ID DHJ-58A-AC23004
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Copper (Cu) | [7440-50-8] | 4,5 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,211 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,606 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,058 | % | | |
| | Lead (Pb) | [7439-92-1] | 0,035 | % | | |
| | Silicon (Si) | [7440-21-3] | 6,58 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,038 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,032 | % | | |
| | Zinc (Zn) | [7440-66-6] | 3,18 | % | | |