

**Aluminium base, RM, Alloy: A380, disc 45x35mm**

Art. ID DHJ-58A-AC23007  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	3,05	%		
	Iron (Fe)	[7439-89-6]	1,68	%		
	Manganese (Mn)	[7439-96-5]	0,189	%		
	Nickel (Ni)	[7440-02-0]	0,189	%		
	Lead (Pb)	[7439-92-1]	0,283	%		
	Silicon (Si)	[7440-21-3]	9,08	%		
	Tin (Sn)	[7440-31-5]	0,299	%		
	Titanium (Ti)	[7440-32-6]	0,193	%		
	Zinc (Zn)	[7440-66-6]	1,09	%		