

Custom Solution - 9 components in Methanol

Art. ID LM24-CUS-70186
Unit 4,5 mL CERTAN® bottle
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Chloropropane	[75-29-6]	1000	mg/L		
	Allyl chloride	[107-05-1]	1000	mg/L		
	Benzene	[71-43-2]	1000	mg/L		
	1,2-Dichloropropane	[78-87-5]	1000	mg/L		
	2,3-Dichloro-1-propene	[78-88-6]	1000	mg/L		
	Epichlorhydrin	[106-89-8]	1000	mg/L		
	cis-1,3-Dichloropropene	[10061-01-5]	1000	mg/L		
	trans-1,3-Dichloropropene	[10061-02-6]	1000	mg/L		
	ne					
	1,2,3-Trichloropropane	[96-18-4]	1000	mg/L		