

Custom Solution in Methanol 4 components

Art. ID LM24-CUS-70194
Unit 2 x 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tetrahydrofuran-d8	[1693-74-9]	200	µg/mL		
	1,4-Dioxane-D8	[17647-74-4]	5000	µg/mL		
	Ethylbenzene-D10	[25837-05-2]	50	µg/mL		
	1,4-Dichlorobenzene D4	[3855-82-1]	100	µg/mL		