

Mix with 5 components in 5% HNO₃

Art. ID LM24-CUS-70196
Unit 2 x 25 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 II E2

Description

Storage conditions: Stored under normal laboratory conditions. Conservation 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Nickel (Ni)	[7440-02-0]	1000	mg/L		
	Copper (Cu)	[7440-50-8]	1000	mg/L		
	Lead (Pb)	[7439-92-1]	1000	mg/L		
	Iron (Fe)	[7439-89-6]	1000	mg/L		
	Zinc (Zn)	[7440-66-6]	1000	mg/L		