

LHKW - Mix 12 in Methanol - Please see LM24-N-10200L-1000ME4x2

Art. ID LM24-N-10240-L-ME4x2
Unit 4 x 2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---|------------|---------------------|-------|--------|--------|
| | Methylene chloride (Dichloromethane) | [75-09-2] | 1000 | µg/mL | | |
| | trans-1,2-Dichloroethene | [156-60-5] | 1000 | µg/mL | | |
| | cis-1,2-Dichloroethene | [156-59-2] | 1000 | µg/mL | | |
| | Chloroform (Trichloromethane) | [67-66-3] | 1000 | µg/mL | | |
| | 1,1,1-Trichloroethane | [71-55-6] | 1000 | µg/mL | | |
| | Carbon tetrachloride (Tetrachloromethane) | [56-23-5] | 1000 | µg/mL | | |
| | Trichloroethylene | [79-01-6] | 1000 | µg/mL | | |
| | Bromodichloromethane | [75-27-4] | 1000 | µg/mL | | |
| | 1,1,2-Trichloroethane | [79-00-5] | 1000 | µg/mL | | |
| | Chlorodibromomethane | [124-48-1] | 1000 | µg/mL | | |
| | Tetrachloroethene | [127-18-4] | 1000 | µg/mL | | |
| | Bromoform (Tribromomethane) | [75-25-2] | 1000 | µg/mL | | |