

**Linolenic acid (C18:3) 1000 µg/mL in Acetone**

Art. ID LM24-N-20000-0115-1000AC5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(Z,Z,Z)-9,12,15-Octadecatrienoic acid (alpha-Linolenic acid) (C18:3 n-3)	[463-40-1]	1000	µg/mL		