

### **Methyl tetracosenoate 100 µg/mL in Methanol**

Art. ID LM24-N-14600-0145-100ME5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	15-Tetracosenoic acid methyl ester (Nervonic acid methyl ester)	[2733-88-2]	100	µg/mL		