

NNK-d4 100 µg/mL in Acetone

Art. ID LM24-N-18000-0788D-100AC5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|---------------|---------------------|-------|--------|--------|
| | NNK-d4 | [764661-24-7] | 100 | µg/mL | | |