

## **Phenanthrene-d10 100 µg/mL in Cyclohexane**

Art. ID LM24-N-14500-1291D-100CY5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1145 II E2

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
|                  | Phenanthrene D10  | [1517-22-2] | 100                 | µg/mL |        |        |