

Phenol-2,3,4,5,6-d5 100 µg/mL in Methanol

Art. ID LM24-N-13400-4345D-100ME5

Unit 5 mL

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------|-------------|---------------------|-------|--------|--------|
| | Phenol-2,3,4,5,6-d5 | [4165-62-2] | 100 | µg/mL | | |