

**Phenol-2,3,4,5,6-d5 100 µg/mL in Methanol**

Art. ID LM24-N-13400-4345D-100ME5

Unit 5 mL

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter   | CAS number  | Concentration/Value | Unit  | Method | Source |
|------------------|---------------------|-------------|---------------------|-------|--------|--------|
|                  | Phenol-2,3,4,5,6-d5 | [4165-62-2] | 100                 | µg/mL |        |        |