

**Tetraethyleneglycol dimethylether 1000 µg/mL in Methanol**

Art. ID LM24-N-13600-0595-1000ME5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter                      | CAS number | Concentration/Value | Unit  | Method | Source |
|------------------|----------------------------------------|------------|---------------------|-------|--------|--------|
|                  | Tetraethylene glycol-di<br>methylether | [143-24-8] | 1000                | µg/mL |        |        |