

Tetrahydrofuran-d8 (THF-d8) 1000 µg/mL in Methanol

Art. ID LM24-N-13600-0601D-1000ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tetrahydrofuran-d8 (THF -d8)	[1693-74-9]	1000	µg/mL		