

Tricosanoic acid (C23:0) 1000 µg/mL in Acetone

Art. ID LM24-N-20000-0140-1000AC5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------|-------------|---------------------|-------|--------|--------|
| | Tricosanoic acid (C23:0) | [2433-96-7] | 1000 | µg/mL | | |