

2,4,7,9-Tetramethyl-5-decyne-4,7-diol (techn.) 1000 µg/mL in Methanol

Art. ID LM24-N-13600-0605-1000ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4,7,9-Tetramethyl-5-decyne-4,7-diol (mix of (±) and meso)	[126-86-3]	1000	µg/mL		