

## **2,4-Dimethyl-3-pentanone 1000 µg/mL in Methanol**

Art. ID LM24-N-13600-0200-1000ME5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4-Dimethyl-3-pentanone	[565-80-0]	1000	µg/mL		
	e					