

## **2,4-Dimethylaniline-d6 100 µg/mL in Methanol**

Art. ID LM24-N-13400-2031D-100ME5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter      | CAS number     | Concentration/Value | Unit  | Method | Source |
|------------------|------------------------|----------------|---------------------|-------|--------|--------|
|                  | 2,4-Dimethylaniline-D6 | [1071170-27-8] | 100                 | µg/mL |        |        |