

2,5-Dichloroaniline-3,4,6-d3 100 µg/mL in Methanol

Art. ID LM24-N-13400-1461D-100ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,5-Dichloroaniline-3,4 ,6-D3	[783321-80-2]	100	µg/mL		