

2,5-Hexanedione-d10 100 µg/mL in Methanol

Art. ID LM24-N-13600-0663D-100ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------|--------------|---------------------|-------|--------|--------|
| | 2,5-Hexanedione-d10 | [97135-07-4] | 100 | µg/mL | | |