

3-Heptanone 1000 µg/mL in Methanol

Art. ID LM24-N-13600-0334-1000ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|-------|--------|--------|
| | 3-Heptanone | [106-35-4] | 1000 | µg/mL | | |