

4,4'-Bis-(dimethylamino)-benzophenone 100 µg/mL in Methanol

Art. ID LM24-N-13400-4320-100ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4,4'-Bis-(dimethylamino)-benzophenone	[90-94-8]	100	µg/mL		