

**4-Toluenesulfonic acid-monohydrate-d7 100 µg/mL in Acetonitrile**

Art. ID LM24-N-15900-4858D-100AN5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Toluenesulfonic acid-monohydrate-d7	[1219795-22-8]	100	µg/mL		