

4-Toluenesulfonic acid-monohydrate-d7 10 µg/mL in Acetone

Art. ID LM24-N-15900-4858D-10AC10
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------------|----------------|---------------------|-------|--------|--------|
| | 4-Toluenesulfonic acid-monohydrate-d7 | [1219795-22-8] | 10 | µg/mL | | |