

**Reagecon Y4 Yellow Colour Reference Solution according to European Pharmacopoeia (EP)**

Art. ID	REA-EPY401
Unit	100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

**Description**

Reagecon's Y4 Yellow Colour Reference Solution is formulated as a reference solution for the examination of the degree of colouration of liquids in accordance with the European Pharmacopoeia Chapter 2.2.2. This product is produced by combining Reagecon's Primary Colour Solutions.