

Density Standard 1,0828g/ml @15°C (ASTM D1480)

| | |
|-----------------|---|
| Art. ID | REA-DEN15160PY |
| Unit | 100 mL |
| Deliverydetails | ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1105 III E1 |

Description

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| at 15°C | Density | | 1,0828 | g/mL | | |