

### **Density Standard 2,0846g/ml @15°C (ASTM D1480)**

Art. ID REA-DEN15240PY  
Unit 100 mL  
Deliverydetails No Dangerous Good /not restricted

#### Description

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| at 15°C          | Density           |            | 2,0846              | g/mL |        |        |