

### **Density Standard 0,7033g/ml @20°C (ASTM D1480)**

Art. ID	REA-DEN20020PY
Unit	100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1262 II E2

#### Description

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 20°C	Density		0,7033	g/mL		