

**Density Standard 1,1651g/ml @20°C (ASTM D1480)**

Art. ID REA-DEN20170PY  
Unit 100 mL  
Deliverydetails No Dangerous Good /not restricted

**Description**

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 20°C	Density		1,1651	g/mL		