

### **Density Standard 0,8084g/ml @25°C (ASTM D1480)**

Art. ID	REA-DEN25090PY
Unit	100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1294 II E2

#### Description

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 25°C	Density		0,8084	g/mL		