

### **Density Standard 1,9108g/ml @25°C (ASTM D1480)**

Art. ID REA-DEN25230PY  
Unit 100 mL  
Deliverydetails No Dangerous Good /not restricted

#### Description

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 25°C	Density		1,9108	g/mL		