

**Aromatic - Mix 19 in Methanol - Please see LM24-N-11750L-100ME2**

Art. ID LM24-N-11760-L-ME2  
Unit 2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene	[71-43-2]	100	µg/mL		
	Toluene	[108-88-3]	100	µg/mL		
	Chlorobenzene	[108-90-7]	100	µg/mL		
	Ethylbenzene	[100-41-4]	100	µg/mL		
	m-Xylene	[108-38-3]	100	µg/mL		
	p-Xylene	[106-42-3]	100	µg/mL		
	o-Xylene	[95-47-6]	100	µg/mL		
	Styrene	[100-42-5]	100	µg/mL		
	Isopropylbenzene (Cumol )	[98-82-8]	100	µg/mL		
	n-Propylbenzene	[103-65-1]	100	µg/mL		
	1,2,3-Trimethylbenzene	[526-73-8]	100	µg/mL		
	1,2,4-Trimethylbenzene	[95-63-6]	100	µg/mL		
	1,3,5-Trimethylbenzene	[108-67-8]	100	µg/mL		
	3-Ethyltoluene	[620-14-4]	100	µg/mL		
	4-Ethyltoluene	[622-96-8]	100	µg/mL		
	2-Ethyltoluene	[611-14-3]	100	µg/mL		
	Indan	[496-11-7]	100	µg/mL		
	Indene	[95-13-6]	100	µg/mL		
	Naphthalene	[91-20-3]	100	µg/mL		