

### **Density Standard 0,8482g/ml @40°C (ASTM D1480)**

Art. ID REA-DEN40110PY  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1294 II E2

#### Description

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 40°C	Density		0,8482	g/mL		