

Density Standard 0,9857g/ml @40°C (ASTM D1480)

Art. ID	REA-DEN40140PY
Unit	100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1105 III E1

Description

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 40°C	Density		0,9857	g/mL		