

### **Density Standard 2,4028g/ml @40°C (ASTM D1480)**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | REA-DEN40260PY                    |
| Unit            | 100 mL                            |
| Deliverydetails | No Dangerous Good /not restricted |

#### Description

Traceable to a National or International standard (eg OIML) // ASTM D1480: Standard Test Method for Density and Relative Density (Specific Gravity) of Viscous Materials by Bingham Pycnometer

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| at 40°C          | Density           |            | 2,4028              | g/mL |        |        |